



# Atomic spectroscopy Workshop

## Interferences

---

Gregory LECORNET  
Atomic spectroscopy product specialist

Agilent Technologies



# Interferences

Interferences are common with any atomic spectroscopy technique, particularly when trying to measure trace-level concentrations in a sample matrix that contains high concentrations of elements that produce line-rich spectra, such as Fe, Ca, and Si.

Sample matrices that are known to generate these types of interferences include geological, metallurgical, and high-matrix environmental samples, such as soils or wastewaters.

The three common types of interferences are

- **physical**
- **chemical**
- **spectral.**

# Physical Interferences

Viscosity

High Dissolved Solids (density)

Acid type or concentration

Surface tension

Organic solvents

# How to check matrix interferences

- Spike recovery (spike volume should be negligible)

$$\%R = 100. \frac{(\text{spiked sample} - \text{sample})}{\text{spike}}$$

- Dilution (>2) recovery

$$\%R = 100. \frac{(\text{dilution factor} \cdot \text{diluted sample})}{\text{sample}}$$

- CRM, certified reference material (matrix should be the same than sample)

# Solving Physical Interference Problems

Dilution (degrades detection limits)

Matrix matching (must be known)

Internal Standardization

Method of Standard Additions

# Internal standard choice

## **IS should:**


- NOT BE present in sample
- NOT BE polluted
- NOT spectrally interfere analyte lines
- NOT be spectrally interfered by sample

For aqueous samples :

For ionic lines: Y, Sc, Ho 1mg/L or In 10 mg/l (in sample mixture)

For in-line addition : Y, Sc, Ho 5mg/L or In 50 mg/l

# Internal Standard Mixing Kit



IS is mixed at a 1 to 5 ratio with Sample



# In-line Internal standard addition (in aqueous)

Total nebulizer flow should remain the same with/without IS  
in-line Addition mixing ratio:

- reduced id ISTD pump tubing to avoid diluting sample
- Ideal relative ratio of ISTD 20-80
- One to one ratio in organics or manual addition

**Mixer:** Y-piece or T-piece

**Clamp pressure:** the clamp pressure is even more important when using ISTD, to be sure ISTD is added with same ratio.





# Chemical Interferences

A sample matrix characteristic which causes an analyte to behave differently in the sample and standard

# Types of Chemical Interferences

Ionization (Na, K, Rb, Cs, Li)

Molecular formation (ie oxides)

Plasma Loading


# Ionization Interferences



As [M] increases, the electron concentration also increases shifting the equilibrium to the left

# EIE


- Concerns easily ionizable elements Li, Na, K, Rb, Cs




Ionization  
Potential  
eV

Li	5,39
Na	5,14
K	4,34
Rb	4,18
Cs	3,89

- ionisation buffer addition 2 g/L  $\text{Cs}(\text{NO}_3)_2$   
(1% for addition in line)




# Ionization Interference Problem #1



For a single element solution of an easily ionized element,  
positive curvature is seen when calibrating

## EIE Study – Effect of high Na on 5 ppm K 766.491nm

	Axial	Axial with Cs	Radial
Na (ppm)	ppm	ppm	ppm
1	5.01	4.91	5.04
10	5.72	4.97	5.03
25	6.61	5.03	4.98
50	7.33	5.02	5.02
100	8.38	5.02	5.03
500	12.18	5.03	5.02
1000	14.48	5.02	4.97




# Solving Chemical Interference Problems

Ionisation buffer : CsCl

Dilution (degrades detection limits)


Matrix Matching

# Use radial Viewing







# Spectral Interferences



# Spectral interference – Cr Interference on Se



-- Se Analyte  
-- Cr interferent

Solution Label	Se 196.026 nm ppm	Se 203.985 nm ppm
Blanc	0.00	0.00
Etalon 1	1.00	1.00
Etalon 2	2.00	2.00
Echantillon 1	1.05	1.17

Solution Label	Se 196.026 nm ppm	Se 203.985 nm ppm
Blanc	0.00	0.00
Etalon 1	1.00	1.00
Etalon 2	2.00	2.00
Echantillon 1	1.05	1.07

Possible interferences on Se (203.985 nm)

Symbol	Wavelength (nm)	Ion	Intensity
Hf	203.909	II	45.4
Re	203.920	I	10.8
Cr	203.932	I	7.2
Ir	203.943	I	35.8
Sb	203.977	I	45.2
Se	203.985	I	161.6
Cr	203.991	II	180.2
Mn	203.998	II	18.6
Ta	204.001	II	45.8
Ta	204.061	II	49.2
Sn	204.066	I	18.4
Fe	204.069	II	43.5